

FUME HOOD BASE CABINETS



Pictured With Optional Recessed Handles

ABOUT

LAS manufactures multiple base cabinets and base stands to support our Fume Hoods and Enclosures, and are designed to accommodate countertops, sinks, plumbing and internal ventilation. Standard base cabinets and stands are available in sitting, standing and ADA compliant heights; made of ½" stress-relieved polypropylene and can be vented for chemical storage. Other non-polypropylene options include flammable storage cabinets and painted metal base stands. All LAS base cabinet structures are manufactured to optimize space and safely support each Fume Hood.



Pictured With Optional Recessed Handles




Pictured With Optional Recessed Handles



Pictured With Optional Recessed Handles

LAS-510V Vented Double Door 3 Foot Fume Hood Base Cabinet With Shelf

(1) Adjustable Shelf with 1" Spill Containment Lip | (1) Polypro Vent 2" FPT Half Coupling



LAS-510V-363528

24" Deep Work Surface Hood

Model

LAS-510V-363522
LAS-510VS-363522

Width Height Depth

36"	35"	22"
36"	35"	22"

30" Deep Work Surface Hood

Model


LAS-510V-363528
LAS-510VS-363528

Width Height Depth

36"	35"	28"
36"	35"	28"

LAS-510V Vented Double Door 4 Foot Fume Hood Base Cabinet With Shelf

(1) Adjustable Shelf with 1" Spill Containment Lip | (1) Polypro Vent 2" FPT Half Coupling



LAS-510VS-483522

24" Deep Work Surface Hood

Model

LAS-510V-483522
LAS-510VS-483522

Width Height Depth

48"	35"	22"
48"	35"	22"

30" Deep Work Surface Hood

Model


LAS-510V-483528
LAS-510VS-483528

Width Height Depth

48"	35"	28"
48"	35"	28"

LAS-525(L/R)V Vented 5 Foot Fume Hood 3 Door Base Cabinet

(1) Divider (Specify Side) | (2) Adjustable Shelves with 1" Spill Containment | (2) Polypro Vents 2" FPT Half Couplings
 **Note: Can Be Ordered With or Without Vents or Shelves



LAS-525LV-603522

24" Deep Work Surface Hood

Model

LAS-520V-603522
LAS-520(L/R)VS-603522

Width Height Depth

60"	35"	22"
60"	35"	22"

30" Deep Work Surface Hood

Model


LAS-510V-603528
LAS-525(L/R)VS-603528

Width Height Depth

60"	35"	28"
60"	35"	28"

LAS-525(L/R)V Vented 6 Foot Fume Hood 3 Door Base Cabinet

(1) Divider (Specify Side) | (2) Adjustable Shelves with 1" Spill Containment | (2) Polypro Vents 2" FPT Half Couplings
 **Note: Can Be Ordered With or Without Vents or Shelves



LAS-525RV-723522


Model

LAS-525V-723522
LAS-525V-723528

Width Height Depth

72"	35"	22"
72"	35"	28"

LAS-530V Vented 8 Foot Fume Hood 4 Door Base Cabinet



(1) Divider | (2) Adjustable Shelves with 1" Spill Containment | (2) Polypro Vents 2" FPT Half Couplings
**Note: Can Be Ordered With or Without Vents or Shelves

Model	Width	Height	Depth
LAS-530V-963522	96"	35"	22"
LAS-530V-963528	96"	35"	28"

LAS-530V-962922