

WALL CABINETS & SHELVING



ABOUT


LAS wall cabinets are single piece construction with single-piece, recessed knuckle hinges and polypropylene pulls. The internal, adjustable shelves are reversible, engineered with a one inch spill containment or a traditional flat surface. All cabinets come with either solid door or polycarbonate glazing, and adjustable shelves with one inch spill containment. Standard wall cabinets are available 12 inches deep and either 24 inches or 30 inches tall. Custom configurations and sizes are available and welcomed.



Pictured With Optional Recessed Handles



LAS-600(L/R) - Hinged Solid Single Door Wall Cabinet




LAS-600L-182412

24" Tall Cabinet comes with (1) Adjustable Shelf with 1" Spill Containment Lip
 30" Tall Cabinet comes with (2) Adjustable Shelves with 1" Spill Containment Lip

Model	Width	Height	Depth	Model	Width	Height	Depth
LAS-600L-182412	18"	24"	12"	LAS-600L-242412	24"	24"	12"
LAS-600R-182412	18"	24"	12"	LAS-600R-242412	24"	24"	12"
LAS-600L-183012	18"	30"	12"	LAS-600L-243012	24"	30"	12"
LAS-600R-183012	18"	30"	12"	LAS-600R-243012	24"	30"	12"

LAS-600(L/R)G Hinged Framed Glazing Single Door Wall Cabinet




LAS-600LG-243012-30

24" Tall Cabinet comes with (1) Adjustable Shelf with 1" Spill Containment Lip
 30" Tall Cabinet comes with (2) Adjustable Shelves with 1" Spill Containment Lip

Model	Width	Height	Depth	Model	Width	Height	Depth
LAS-600LG-182412	18"	24"	12"	LAS-600LG-242412	24"	24"	12"
LAS-600RG-182412	18"	24"	12"	LAS-600RG-242412	24"	24"	12"
LAS-600LG-183012	18"	30"	12"	LAS-600LG-243012	24"	30"	12"
LAS-600RG-183012	18"	30"	12"	LAS-600RG-243012	24"	30"	12"

LAS-610 Hinged Solid Double Door Wall Cabinet




LAS-610-302412

24" Tall Cabinet comes with (1) Adjustable Shelf with 1" Spill Containment Lip
 30" Tall Cabinet comes with (2) Adjustable Shelves with 1" Spill Containment Lip

Model	Width	Height	Depth	Model	Width	Height	Depth
LAS-610-302412	30"	24"	12"	LAS-610-303012	30"	30"	12"
LAS-610-362412	36"	24"	12"	LAS-610-363012	36"	30"	12"
LAS-610-422412	42"	24"	12"	LAS-610-423012	42"	30"	12"
LAS-610-482412	48"	24"	12"	LAS-610-483012	48"	30"	12"

LAS-620 Solid Sliding Double Door Wall Cabinet




LAS-620-303012

24" Tall Cabinet comes with (1) Adjustable Shelf with 1" Spill Containment Lip
 30" Tall Cabinet comes with (2) Adjustable Shelves with 1" Spill Containment Lip
 Note: Recessed handles not available on sliding doors

Model	Width	Height	Depth	Model	Width	Height	Depth
LAS-620-302412	30"	24"	12"	LAS-620-303012	30"	30"	12"
LAS-620-362412	36"	24"	12"	LAS-620-363012	36"	30"	12"
LAS-620-422412	42"	24"	12"	LAS-620-423012	42"	30"	12"
LAS-620-482412	48"	24"	12"	LAS-620-483012	48"	30"	12"

LAS-610G Hinged Framed Glazing Double Door Wall Cabinet




24" Tall Cabinet comes with (1) Adjustable Shelf with 1" Spill Containment Lip
 30" Tall Cabinet comes with (2) Adjustable Shelves with 1" Spill Containment Lip

Model	Width	Height	Depth	Model	Width	Height	Depth
LAS-610G-302412	30"	24"	12"	LAS-610G-303012	30"	30"	12"
LAS-610G-362412	36"	24"	12"	LAS-610G-363012	36"	30"	12"
LAS-610G-422412	42"	24"	12"	LAS-610G-423012	42"	30"	12"
LAS-610G-482412	48"	24"	12"	LAS-610G-483012	48"	30"	12"

LAS-610G-303012

LAS-620G Framed Glazing Sliding Double Door Wall Cabinet




24" Tall Cabinet comes with (1) Adjustable Shelf with 1" Spill Containment Lip
 30" Tall Cabinet comes with (2) Adjustable Shelves with 1" Spill Containment Lip
 Note: Recessed handles not available on sliding doors

Model	Width	Height	Depth	Model	Width	Height	Depth
LAS-620G-302412	30"	24"	12"	LAS-620G-303012	30"	30"	12"
LAS-620G-362412	36"	24"	12"	LAS-620G-363012	36"	30"	12"
LAS-620G-422412	42"	24"	12"	LAS-620G-423012	42"	30"	12"
LAS-620G-482412	48"	24"	12"	LAS-620G-483012	48"	30"	12"

LAS-620G-363012

LAS-610(L/R)B Hinged Solid Double Door Blind Corner Wall Cabinet




24" Tall Cabinet comes with (1) Adjustable Shelf with 1" Spill Containment Lip
 30" Tall Cabinet comes with (2) Adjustable Shelves with 1" Spill Containment Lip
 Note: Specify Left or Right Blind Corner area.

Model	Width	Height	Depth	Model	Width	Height	Depth
LAS-610LB-302412	42"	24"	12"	LAS-610LB-303012	42"	30"	12"
LAS-610RB-302412	42"	24"	12"	LAS-610RB-303012	42"	30"	12"
LAS-610LB-362412	48"	24"	12"	LAS-610LB-363012	48"	30"	12"
LAS-610RB-362412	48"	24"	12"	LAS-610RB-363012	48"	30"	12"
LAS-610LB-422412	54"	24"	12"	LAS-610LB-423012	54"	30"	12"
LAS-610RB-422412	54"	24"	12"	LAS-610RB-423012	54"	30"	12"
LAS-610LB-482412	60"	24"	12"	LAS-610LB-483012	60"	30"	12"
LAS-610RB-482412	60"	24"	12"	LAS-610RB-483012	60"	30"	12"

LAS-610LB-302412

LAS-610(L/R)BG Hinged Framed Glazing Double Door Blind Corner Wall Cabinet


24" Tall Cabinet comes with (1) Fixed Shelf with 1" Spill Containment Lip
 30" Tall Cabinet comes with (2) Fixed Shelves with 1" Spill Containment Lip
 Note: Specify Left or Right Blind Corner area.



LAS-610LBG-303012

Model	Width	Height	Depth	Model	Width	Height	Depth
LAS-610LBG-302412	42"	24"	12"	LAS-610LBG-303012	42"	30"	12"
LAS-610RBG-302412	42"	24"	12"	LAS-610RBG-303012	42"	30"	12"
LAS-610LBG-362412	48"	24"	12"	LAS-610LBG-363012	48"	30"	12"
LAS-610RBG-362412	48"	24"	12"	LAS-610RBG-363012	48"	30"	12"
LAS-610LBG-422412	54"	24"	12"	LAS-610LBG-423012	54"	30"	12"
LAS-610RBG-422412	54"	24"	12"	LAS-610RBG-423012	54"	30"	12"
LAS-610LBG-482412	60"	24"	12"	LAS-610LBG-483012	60"	30"	12"
LAS-610RBG-482412	60"	24"	12"	LAS-610RBG-483012	60"	30"	12"

LAS-610N Open Shelf Wall Cabinet




LAS-610N-18 & LAS-610N-24 Models: (2) Adjustable Shelves with 1" Spill Containment Lip
 All Other Models: (4) Adjustable Shelves with 1" Spill Containment Lip

LAS-610N-422412

Model	Width	Height	Depth	Model	Width	Height	Depth
LAS-610N-182412	18"	24"	12"	LAS-610N-183012	18"	30"	12"
LAS-610N-242412	24"	24"	12"	LAS-610N-243012	24"	30"	12"
LAS-610N-302412	30"	24"	12"	LAS-610N-303012	30"	30"	12"
LAS-610N-362412	36"	24"	12"	LAS-610N-363012	36"	30"	12"
LAS-610N-422412	42"	24"	12"	LAS-610N-423012	42"	30"	12"
LAS-610N-482412	48"	24"	12"	LAS-610N-483012	48"	30"	12"

LAS-640N Single Open Shelf Corner Wall Cabinet



24" Tall Cabinet comes with (1) Fixed Shelf with 1" Spill Containment Lip
 30" Tall Cabinet comes with (2) Fixed Shelves with 1" Spill Containment Lip

LAS-640N-163012

Model	Width	Height	Depth	Model	Width	Height	Depth
LAS-640N-122412	12"	24"	12"	LAS-640N-123012	12"	30"	12"
LAS-640N-162412	16"	24"	12"	LAS-640N-163012	16"	30"	12"



9650 Newton Ave. S
Bloomington, MN 55431



877.410.4910



info@labairesystems.com



www.labairesystems.com