

## BASE CABINETS WITH DRAWERS & DOORS




Pictured With Optional Recessed Handles

### ABOUT

LAS drawer/door combination base cabinets come with multiple standard configurations of doors(single or double) with up to four drawers. All cabinet doors are fastened with single-piece, recessed knuckle hinges and feature single piece polypropylene pulls. The internal shelves are removeable, adjustable by height and reversible, engineered with a one inch spill containment or a traditional flat surface. Standard base cabinets are available in sitting, standing and ADA compliant heights. LAS uses low friction drawer guides/glides material for effortless operation. LAS also offers an optional soft-close mechanism, allowing the drawers to close silently and smoothly. Custom configurations and sizes are available and welcomed.



## LAS-400(L/R) Hinged Solid Single Left/Right — Door/Drawer Base Cabinet



LAS-400R-183522

(1) Adjustable Shelf with 1" Spill Containment Lip

### Sitting Height

Model
LAS-400L-182922
LAS-400R-182922
LAS-400L-242922
LAS-400R-242922


Width	Height	Depth
18"	29"	22"
18"	29"	22"
24"	29"	22"
24"	29"	22"

### Standing Height

Model
LAS-400L-183522
LAS-400R-183522
LAS-400L-243522
LAS-400R-243522

Width	Height	Depth
18"	35"	22"
18"	35"	22"
24"	35"	22"
24"	35"	22"

## LAS-405 Hinged Solid Double Door With Single Drawer Base Cabinet



LAS-405-363522

(1) Adjustable Shelf with 1" Spill Containment Lip

### Sitting Height

Model
LAS-405-302922
LAS-405-362922
LAS-405-422922
LAS-405-482922


Width	Height	Depth
30"	29"	22"
36"	29"	22"
42"	29"	22"
48"	29"	22"

### Standing Height

Model
LAS-405-303522
LAS-405-363522
LAS-405-423522
LAS-405-483522

Width	Height	Depth
30"	35"	22"
36"	35"	22"
42"	35"	22"
48"	35"	22"

## LAS-410 Hinged Solid Double Door With (2) Drawer Base Cabinet



LAS-410-483522

(1) Adjustable Shelf with 1" Spill Containment Lip

### Sitting Height

Model
LAS-410-302922
LAS-410-362922
LAS-410-422922
LAS-410-482922

Width	Height	Depth
30"	29"	22"
36"	29"	22"
42"	29"	22"
48"	29"	22"


### Standing Height

Model
LAS-410-303522
LAS-410-363522
LAS-410-423522
LAS-410-483522

Width	Height	Depth
30"	35"	22"
36"	35"	22"
42"	35"	22"
48"	35"	22"

## LAS-401(L/R) Hinged Single Left/Right Door/ (3) Drawer Base Cabinet

(1) Adjustable Shelf with 1" Spill Containment Lip



LAS-401L-483522

### Sitting Height


Model	Width	Height	Depth
LAS-401L-302922	30"	29"	22"
LAS-401R-302922	30"	29"	22"
LAS-401L-362922	36"	29"	22"
LAS-401R-362922	36"	29"	22"
LAS-401L-422922	42"	29"	22"
LAS-401R-422922	42"	29"	22"
LAS-401L-482922	48"	29"	22"
LAS-401R-482922	48"	29"	22"

### Standing Height

Model	Width	Height	Depth
LAS-401L-303522	30"	35"	22"
LAS-401R-303522	30"	35"	22"
LAS-401L-363522	36"	35"	22"
LAS-401R-363522	36"	35"	22"
LAS-401L-423522	42"	35"	22"
LAS-401R-423522	42"	35"	22"
LAS-401L-483522	48"	35"	22"
LAS-401R-483522	48"	35"	22"

## LAS-402(L/R) Hinged Single Left/Right Door/ (4) Drawer Base Cabinet

(1) Adjustable Shelf with 1" Spill Containment Lip



LAS-402R-483522

### Sitting Height

Model	Width	Height	Depth
LAS-402L-302922	30"	29"	22"
LAS-402R-302922	30"	29"	22"
LAS-402L-362922	36"	29"	22"
LAS-402R-362922	36"	29"	22"
LAS-402L-422922	42"	29"	22"
LAS-402R-422922	42"	29"	22"
LAS-402L-482922	48"	29"	22"
LAS-402R-482922	48"	29"	22"

### Standing Height

Model	Width	Height	Depth
LAS-402L-303522	30"	35"	22"
LAS-402R-303522	30"	35"	22"
LAS-402L-363522	36"	35"	22"
LAS-402R-363522	36"	35"	22"
LAS-402L-423522	42"	35"	22"
LAS-402R-423522	42"	35"	22"
LAS-402L-483522	48"	35"	22"
LAS-402R-483522	48"	35"	22"