

BASE, SINK, VENTED, AND TALL CABINETS



ABOUT

LAS manufactures multiple variations of standard base cabinets to accommodate countertops, sinks, plumbing, internal ventilation, and fume hoods. All cabinet doors are fastened with single-piece, recessed knuckle hinges and single piece polypropylene pulls. The internal, adjustable shelves are reversible, engineered for a one inch spill containment or a traditional flat surface. Standard base cabinets are available in sitting, standing and ADA compliant heights. The LAS tall cabinets range in height from five to seven feet high. Custom configurations and sizes are available and welcomed.



Pictured With Optional Recessed Handles




Pictured With Optional Recessed Handles



Pictured With Optional Recessed Handles

LAS-500(L/R) Left/Right Hinged Full Door Base Cabinet With Shelf

Left or Right Hinged | (1) Adjustable Shelf with 1" Spill Containment Lip



LAS-500R-182922

Sitting Height

Model
LAS-500L -182922
LAS-500R-182922
LAS-500L-242922
LAS-500R-242922

Width	Height	Depth
18"	29"	22"
18"	29"	22"
24"	29"	22"
24"	29"	22"


Standing Height

Model
LAS-500L -183522
LAS-500R-183522
LAS-500L-243522
LAS-500R-243522

Width	Height	Depth
18"	35"	22"
18"	35"	22"
24"	35"	22"
24"	35"	22"

LAS-510 Hinged Double Door Base Cabinet With Shelf

(1) Adjustable Shelf with 1" Spill Containment Lip



LAS-510-302922

Sitting Height

Model
LAS-510-302922
LAS-510-362922
LAS-510-422922
LAS-510-482922


Width	Height	Depth
30"	29"	22"
36"	29"	22"
42"	29"	22"
48"	29"	22"

Standing Height

Model
LAS-510-303522
LAS-510-363522
LAS-510-423522
LAS-510-483522

Width	Height	Depth
30"	35"	22"
36"	35"	22"
42"	35"	22"
48"	35"	22"

LAS-510S Hinged Double Door Sink Cabinet No Shelf



LAS-510S-482922

Sitting Height

Model
LAS-510S-302922
LAS-510S-362922
LAS-510S-422922
LAS-510S-482922

Width	Height	Depth
30"	29"	22"
36"	29"	22"
42"	29"	22"
48"	29"	22"


Standing Height

Model
LAS-510S-303522
LAS-510S-363522
LAS-510S-423522
LAS-510S-483522

Width	Height	Depth
30"	35"	22"
36"	35"	22"
42"	35"	22"
48"	35"	22"

LAS-510V Vented Acid Storage Double Door Base Cabinet With Shelf

(1) Adjustable Shelf with 1" Spill Containment Lip | (1) Polypro Vent 2" FPT Half Coupling



Sitting Height

Model
LAS-510V-302922
LAS-510V-362922
LAS-510V-422922
LAS-510V-482922


Width	Height	Depth
30"	29"	22"
36"	29"	22"
42"	29"	22"
48"	29"	22"

Standing Height

Model
LAS-510V-303522
LAS-510V-363522
LAS-510V-423522
LAS-510V-483522

Width	Height	Depth
30"	35"	22"
36"	35"	22"
42"	35"	22"
48"	35"	22"

LAS-511S ADA Open Sink Cabinet With Angle Board



LAS-511S-483322


- (1) Removable Angle Board
- (1) Removable Plumbing Access Panel

Model

- LAS-511S-303322
- LAS-511S-363322
- LAS-511S-423322
- LAS-511S-483322

Width	Height	Depth
30"	33"	22"
36"	33"	22"
42"	33"	22"
48"	33"	22"

LAS-500(L/R) Left/Right Solid Single Door 5 Foot Tall Cabinet



LAS-500L-246018


- Left or Right Hinged
- Shelves Have 1" Spill Containment Lip
- (1) Fixed Shelf | (2) Adjustable Shelves

Model

- LAS-500L -186018
- LAS-500R-186018
- LAS-500L-246018
- LAS-500R-246018

Width	Height	Depth
18"	60"	18"
18"	60"	18"
24"	60"	18"
24"	60"	18"

LAS-500(L/R) Left/Right Solid Single Door 7 Foot Tall Cabinet



LAS-500L-248422


- Left or Right Hinged
- Shelves Have 1" Spill Containment Lip
- (1) Fixed Shelf | (4) Adjustable Shelves

Model

- LAS-500L -188422
- LAS-500R-188422
- LAS-500L-248422
- LAS-500R-248422

Width	Height	Depth
18"	84"	22"
18"	84"	22"
24"	84"	22"
24"	84"	22"

LAS-500(L/R)G Left/Right Framed Glazing Single Door 5 Foot Tall Cabinet



LAS-500RG-246018


- Left or Right Hinged
- Shelves Have 1" Spill Containment Lip
- (1) Fixed Shelf | (2) Adjustable Shelves

Model

- LAS-500LG -186018
- LAS-500RG-186018
- LAS-500LG-246018
- LAS-500RG-246018

Width	Height	Depth
18"	60"	18"
18"	60"	18"
24"	60"	18"
24"	60"	18"

LAS-500(L/R)G Left/Right Framed Glazing Single Door 7 Foot Tall Cabinet




Left or Right Hinged
Shelves Have 1" Spill Containment Lip
(1) Fixed Shelf | (4) Adjustable Shelves

Model	Width	Height	Depth
LAS-500LG-188422	18"	84"	22"
LAS-500RG-188422	18"	84"	22"
LAS-500LG-248422	24"	84"	22"
LAS-500RG-248422	24"	84"	22"

LAS-500RG-188422

LAS-500(L/R)V Vented Left/Right Single Door 5 Foot Tall Cabinet




Left or Right Hinged
Shelves Have 1" Spill Containment Lip
(1) Fixed Shelf | (2) Adjustable Shelves | (1) 2" FPT Vent Half Coupling

Model	Width	Height	Depth
LAS-500LV-186018	18"	60"	18"
LAS-500RV-186018	18"	60"	18"
LAS-500LV-246018	24"	60"	18"
LAS-500RV-246018	24"	60"	18"

LAS-500LV-246018

LAS-500(L/R)V Vented Left/Right Single Door 7 Foot Tall Cabinet




Left or Right Hinged
Shelves Have 1" Spill Containment Lip
(1) Fixed Shelf | (4) Adjustable Shelves | (1) 2" FPT Vent Half Coupling

Model	Width	Height	Depth
LAS-500LV-188422	18"	84"	22"
LAS-500RV-188422	18"	84"	22"
LAS-500LV-248422	24"	84"	22"
LAS-500RV-248422	24"	84"	22"

LAS-500LV-248422

LAS-510 Hinged Solid Double Door 5 Foot Tall Cabinet

Shelves Have 1" Spill Containment Lip | (1) Fixed Shelf | (2) Adjustable Shelves




Model	Width	Height	Depth
LAS-510-306018	30"	60"	18"
LAS-510-366018	36"	60"	18"
LAS-510-426018	42"	60"	18"
LAS-510-486018	48"	60"	18"

LAS-510-486018

LAS-510 Hinged Solid Double Door 7 Foot Tall Cabinet

Shelves Have 1" Spill Containment Lip | (1) Fixed Shelf | (4) Adjustable Shelves



LAS-510-488422


Model

LAS-510-308422
 LAS-510-368422
 LAS-510-428422
 LAS-510-488422

Width	Height	Depth
30"	84"	22"
36"	84"	22"
42"	84"	22"
48"	84"	22"

LAS-510G Hinged Framed Glazing Double Door 5 Foot Tall Cabinet

Shelves Have 1" Spill Containment Lip | (1) Fixed Shelf | (2) Adjustable Shelves



LAS-510-488418


Model

LAS-510G-306018
 LAS-510G-366018
 LAS-510G-426018
 LAS-510G-486018

Width	Height	Depth
30"	60"	18"
36"	60"	18"
42"	60"	18"
48"	60"	18"

LAS-510G Hinged Framed Glazing Double Door 7 Foot Tall Cabinet

Shelves Have 1" Spill Containment Lip | (1) Fixed Shelf | (4) Adjustable Shelves



LAS-510G-308422


Model

LAS-510G-308422
 LAS-510G-368422
 LAS-510G-428422
 LAS-510G-488422

Width	Height	Depth
30"	84"	22"
36"	84"	22"
42"	84"	22"
48"	84"	22"

LAS-510V Hinged Vented Double Door 5 Foot Tall Cabinet

Shelves Have 1" Spill Containment Lip | (1) Fixed Shelf | (2) Adjustable Shelves



LAS-510V-426018


Model

LAS-510V-306018
 LAS-510V-366018
 LAS-510V-426018
 LAS-510V-486018

Width	Height	Depth
30"	60"	18"
36"	60"	18"
42"	60"	18"
48"	60"	18"

LAS-510V Hinged Vented Double Door 7 Foot Tall Cabinet

Shelves Have 1" Spill Containment Lip | (1) Fixed Shelf | (4) Adjustable Shelves



LAS-510V-428422


Model

LAS-510V-308422
 LAS-510V-368422
 LAS-510V-428422
 LAS-510V-488422

Width	Height	Depth
30"	84"	22"
36"	84"	22"
42"	84"	22"
48"	84"	22"

LAS-510N Open Shelf Bookcase Base Cabinet

LAS-510N-18 & LAS-510N-24 Models: (2) Adjustable Shelves with 1" Spill Containment Lip
 All Other Models: (4) Adjustable Shelves with 1" Spill Containment Lip



LAS-510N-423522

Sitting Height

Model

LAS-510N-182922
 LAS-510N-242922
 LAS-510N-302922
 LAS-510N-362922
 LAS-510N-422922
 LAS-510N-482922

Width	Height	Depth
18"	29"	22"
24"	29"	22"
30"	29"	22"
36"	29"	22"
42"	29"	22"
48"	29"	22"

Standing Height


Model

LAS-510N-183522
 LAS-510N-243522
 LAS-510N-303522
 LAS-510N-363522
 LAS-510N-423522
 LAS-510N-483522

Width	Height	Depth
18"	35"	22"
24"	35"	22"
30"	35"	22"
36"	35"	22"
42"	35"	22"
48"	35"	22"

LAS-510(L/R)B Hinged Solid Double Door Left/Right Blind Corner Cabinet

(1) Fixed Shelf with 1" spill Containment Lip



LAS-510LB-363522

Sitting Height

Model

LAS-510LB-302922
 LAS-510RB-302922
 LAS-510LB-362922
 LAS-510RB-362922
 LAS-510LB-422922
 LAS-510RB-422922
 LAS-510LB-482922
 LAS-510RB-482922

Width	Height	Depth
52"	29"	22"
52"	29"	22"
58"	29"	22"
58"	29"	22"
64"	29"	22"
64"	29"	22"
70"	29"	22"
70"	29"	22"


Standing Height

Model

LAS-510LB-303522
 LAS-510RB-303522
 LAS-510LB-363522
 LAS-510RB-363522
 LAS-510LB-423522
 LAS-510EB-423522
 LAS-510LB-483522
 LAS-510RB-483522

Width	Height	Depth
52"	35"	22"
52"	35"	22"
58"	35"	22"
58"	35"	22"
64"	35"	22"
64"	35"	22"
70"	35"	22"
70"	35"	22"

LAS-540(L/R) Hinged Solid Single Left or Right Door Corner Base Cabinet



LAS-540R-162922

Sitting Height

Model

LAS-540L-122922
 LAS-540R-122922
 LAS-540L-162922
 LAS-540R-162922

Width	Height	Depth
12"	29"	22"
12"	29"	22"
16"	29"	22"
16"	29"	22"

Standing Height

Model

LAS-540L-123522
 LAS-540R-123522
 LAS-540L-163522
 LAS-540R-163522

Width	Height	Depth
12"	35"	22"
12"	35"	22"
16"	35"	22"
16"	35"	22"