

FUME HOOD BASE CABINETS



Pictured With Optional Recessed Handles

ABOUT

LAS manufactures multiple base cabinets and base stands to support our Fume Hoods and Enclosures, and are designed to accommodate countertops, sinks, plumbing and internal ventilation. Standard base cabinets and stands are available in sitting, standing and ADA compliant heights; made of ½" stress-relieved polypropylene and can be vented for chemical storage. Other non-polypropylene options include flammable storage cabinets and painted metal base stands. All LAS base cabinet structures are manufactured to optimize space and safely support each Fume Hood.



Pictured With Optional Recessed Handles




Pictured With Optional Recessed Handles



Pictured With Optional Recessed Handles

LAS-510V Vented Double Door 3 Foot Fume Hood Base Cabinet With Shelf

(1) Adjustable Shelf with 1" Spill Containment Lip | (1) Polypro Vent 2" FPT Half Coupling



LAS-510V-363528

24" Deep Work Surface Hood

Model

LAS-510V-363522

LAS-510VS-363522

Width	Height	Depth
36"	35"	22"
36"	35"	22"

30" Deep Work Surface Hood

Model


LAS-510V-363528

LAS-510VS-363528

Width	Height	Depth
36"	35"	28"
36"	35"	28"

LAS-510V Vented Double Door 4 Foot Fume Hood Base Cabinet With Shelf

(1) Adjustable Shelf with 1" Spill Containment Lip | (1) Polypro Vent 2" FPT Half Coupling



LAS-510VS-483522

24" Deep Work Surface Hood

Model

LAS-510V-483522

LAS-510VS-483522

Width	Height	Depth
48"	35"	22"
48"	35"	22"

30" Deep Work Surface Hood

Model

LAS-510V-483528


LAS-510VS-483528

Width	Height	Depth
48"	35"	28"
48"	35"	28"

LAS-525(L/R)V Vented 5 Foot Fume Hood 3 Door Base Cabinet

(1) Divider (Specify Side) | (2) Adjustable Shelves with 1" Spill Containment | (2) Polypro Vents 2" FPT Half Couplings

**Note: Can Be Ordered With or Without Vents or Shelves



LAS-525LV-603522

24" Deep Work Surface Hood

Model

LAS-520V-603522

LAS-520(L/R)VS-603522

Width	Height	Depth
60"	35"	22"
60"	35"	22"

30" Deep Work Surface Hood

Model

LAS-510V-603528


LAS-525(L/R)VS-603528

Width	Height	Depth
60"	35"	28"
60"	35"	28"

LAS-525(L/R)V Vented 6 Foot Fume Hood 3 Door Base Cabinet

(1) Divider (Specify Side) | (2) Adjustable Shelves with 1" Spill Containment | (2) Polypro Vents 2" FPT Half Couplings

**Note: Can Be Ordered With or Without Vents or Shelves



LAS-525RV-723522

Model

LAS-525V-723522


LAS-525V-723528

Width	Height	Depth
72"	35"	22"
72"	35"	28"

LAS-530V Vented 8 Foot Fume Hood 4 Door Base Cabinet

(1) Divider | (2) Adjustable Shelves with 1" Spill Containment | (2) Polypro Vents 2" FPT Half Couplings

**Note: Can Be Ordered With or Without Vents or Shelves



LAS-530V-962922

Model

LAS-530V-963522

LAS-530V-963528

Width	Height	Depth
96"	35"	22"
96"	35"	28"