

## WALL CABINETS & SHELVING




### ABOUT

LAS wall cabinets are single piece construction with single-piece, recessed knuckle hinges and polypropylene pulls. The internal, adjustable shelves are reversible, engineered with a one inch spill containment or a traditional flat surface. All cabinets come with either solid door or polycarbonate glazing, and adjustable shelves with one inch spill containment. Standard wall cabinets are available 12 inches deep and either 24 inches or 30 inches tall. Custom configurations and sizes are available and welcomed.



## LAS-600(L/R) - Hinged Solid Single Door Wall Cabinet



LAS-600L-182412

24" Tall Cabinet comes with (1) Adjustable Shelf with 1" Spill Containment Lip  
 30" Tall Cabinet comes with (2) Adjustable Shelves with 1" Spill Containment Lip


Model
LAS-600L-182412
LAS-600R-182412
LAS-600L-183012
LAS-600R-183012

Width	Height	Depth
18"	24"	12"
18"	24"	12"
18"	30"	12"
18"	30"	12"

Model
LAS-600L-242412
LAS-600R-242412
LAS-600L-243012
LAS-600R-243012

Width	Height	Depth
24"	24"	12"
24"	24"	12"
24"	30"	12"
24"	30"	12"

## LAS-600(L/R)G Hinged Framed Glazing Single Door Wall Cabinet



LAS-600LG-243012-30

24" Tall Cabinet comes with (1) Adjustable Shelf with 1" Spill Containment Lip  
 30" Tall Cabinet comes with (2) Adjustable Shelves with 1" Spill Containment Lip


Model
LAS-600LG-182412
LAS-600RG-182412
LAS-600LG-183012
LAS-600RG-183012

Width	Height	Depth
18"	24"	12"
18"	24"	12"
18"	30"	12"
18"	30"	12"

Model
LAS-600LG-242412
LAS-600RG-242412
LAS-600LG-243012
LAS-600RG-243012

Width	Height	Depth
24"	24"	12"
24"	24"	12"
24"	30"	12"
24"	30"	12"

## LAS-610 Hinged Solid Double Door Wall Cabinet



LAS-610-302412

24" Tall Cabinet comes with (1) Adjustable Shelf with 1" Spill Containment Lip  
 30" Tall Cabinet comes with (2) Adjustable Shelves with 1" Spill Containment Lip


Model
LAS-610-302412
LAS-610-362412
LAS-610-422412
LAS-610-482412

Width	Height	Depth
30"	24"	12"
36"	24"	12"
42"	24"	12"
48"	24"	12"

Model
LAS-610-303012
LAS-610-363012
LAS-610-423012
LAS-610-483012

Width	Height	Depth
30"	30"	12"
36"	30"	12"
42"	30"	12"
48"	30"	12"

## LAS-620 Solid Sliding Double Door Wall Cabinet



LAS-620-303012

24" Tall Cabinet comes with (1) Adjustable Shelf with 1" Spill Containment Lip  
 30" Tall Cabinet comes with (2) Adjustable Shelves with 1" Spill Containment Lip  
 Note: Recessed handles not available on sliding doors


Model
LAS-620-302412
LAS-620-362412
LAS-620-422412
LAS-620-482412

Width	Height	Depth
30"	24"	12"
36"	24"	12"
42"	24"	12"
48"	24"	12"

Model
LAS-620-303012
LAS-620-363012
LAS-620-423012
LAS-620-483012

Width	Height	Depth
30"	30"	12"
36"	30"	12"
42"	30"	12"
48"	30"	12"

## LAS-610G Hinged Framed Glazing Double Door Wall Cabinet



24" Tall Cabinet comes with (1) Adjustable Shelf with 1" Spill Containment Lip  
 30" Tall Cabinet comes with (2) Adjustable Shelves with 1" Spill Containment Lip

Model
LAS-610G-302412
LAS-610G-362412
LAS-610G-422412
LAS-610G-482412


Width	Height	Depth
30"	24"	12"
36"	24"	12"
42"	24"	12"
48"	24"	12"

Model
LAS-610G-303012
LAS-610G-363012
LAS-610G-423012
LAS-610G-483012

Width	Height	Depth
30"	30"	12"
36"	30"	12"
42"	30"	12"
48"	30"	12"

LAS-610G-303012

## LAS-620G Framed Glazing Sliding Double Door Wall Cabinet



24" Tall Cabinet comes with (1) Adjustable Shelf with 1" Spill Containment Lip  
 30" Tall Cabinet comes with (2) Adjustable Shelves with 1" Spill Containment Lip  
 Note: Recessed handles not available on sliding doors

Model
LAS-620G-302412
LAS-620G-362412
LAS-620G-422412
LAS-620G-482412


Width	Height	Depth
30"	24"	12"
36"	24"	12"
42"	24"	12"
48"	24"	12"

Model
LAS-620G-303012
LAS-620G-363012
LAS-620G-423012
LAS-620G-483012

Width	Height	Depth
30"	30"	12"
36"	30"	12"
42"	30"	12"
48"	30"	12"

LAS-620G-363012

## LAS-610(L/R)B Hinged Solid Double Door Blind Corner Wall Cabinet



24" Tall Cabinet comes with (1) Adjustable Shelf with 1" Spill Containment Lip  
 30" Tall Cabinet comes with (2) Adjustable Shelves with 1" Spill Containment Lip  
 Note: Specify Left or Right Blind Corner area.

Model
LAS-610LB-302412
LAS-610RB-302412
LAS-610LB-362412
LAS-610RB-362412
LAS-610LB-422412
LAS-610RB-422412
LAS-610LB-482412
LAS-610RB-482412

Width	Height	Depth
42"	24"	12"
42"	24"	12"
48"	24"	12"
48"	24"	12"
54"	24"	12"
54"	24"	12"
60"	24"	12"
60"	24"	12"


Model
LAS-610LB-303012
LAS-610RB-303012
LAS-610LB-363012
LAS-610RB-363012
LAS-610LB-423012
LAS-610RB-423012
LAS-610LB-483012
LAS-610RB-483012

Width	Height	Depth
42"	30"	12"
42"	30"	12"
48"	30"	12"
48"	30"	12"
54"	30"	12"
54"	30"	12"
60"	30"	12"
60"	30"	12"

LAS-610LB-302412

## LAS-610(L/R)BG Hinged Framed Glazing Double Door Blind Corner Wall Cabinet

24" Tall Cabinet comes with (1) Fixed Shelf with 1" Spill Containment Lip  
 30" Tall Cabinet comes with (2) Fixed Shelves with 1" Spill Containment Lip  
 Note: Specify Left or Right Blind Corner area.




LAS-610L BG-303012

Model	Width	Height	Depth	Model	Width	Height	Depth
LAS-610L BG-302412	42"	24"	12"	LAS-610L BG-303012	42"	30"	12"
LAS-610R BG-302412	42"	24"	12"	LAS-610R BG-303012	42"	30"	12"
LAS-610L BG-362412	48"	24"	12"	LAS-610L BG-363012	48"	30"	12"
LAS-610R BG-362412	48"	24"	12"	LAS-610R BG-363012	48"	30"	12"
LAS-610L BG-422412	54"	24"	12"	LAS-610L BG-423012	54"	30"	12"
LAS-610R BG-422412	54"	24"	12"	LAS-610R BG-423012	54"	30"	12"
LAS-610L BG-482412	60"	24"	12"	LAS-610L BG-483012	60"	30"	12"
LAS-610R BG-482412	60"	24"	12"	LAS-610R BG-483012	60"	30"	12"

## LAS-610N Open Shelf Wall Cabinet

LAS-610N-18 & LAS-610N-24 Models: (2) Adjustable Shelves with 1" Spill Containment Lip  
 All Other Models: (4) Adjustable Shelves with 1" Spill Containment Lip




LAS-610N-422412

Model	Width	Height	Depth	Model	Width	Height	Depth
LAS-610N-182412	18"	24"	12"	LAS-610N-183012	18"	30"	12"
LAS-610N-242412	24"	24"	12"	LAS-610N-243012	24"	30"	12"
LAS-610N-302412	30"	24"	12"	LAS-610N-303012	30"	30"	12"
LAS-610N-362412	36"	24"	12"	LAS-610N-363012	36"	30"	12"
LAS-610N-422412	42"	24"	12"	LAS-610N-423012	42"	30"	12"
LAS-610N-482412	48"	24"	12"	LAS-610N-483012	48"	30"	12"

## LAS-640N Single Open Shelf Corner Wall Cabinet

24" Tall Cabinet comes with (1) Fixed Shelf with 1" Spill Containment Lip  
 30" Tall Cabinet comes with (2) Fixed Shelves with 1" Spill Containment Lip



LAS-640N-163012

Model	Width	Height	Depth	Model	Width	Height	Depth
LAS-640N-122412	12"	24"	12"	LAS-640N-123012	12"	30"	12"
LAS-640N-162412	16"	24"	12"	LAS-640N-163012	16"	30"	12"

DISTRIBUTED BY:

